混凝试验搅拌机
Coagulation test mixer

ZKS- I -6L

产品使用说明
Product User Manual

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CONTENTS

Preface	1
Product Overview	1
Product Introduction	1
Main role	1
Scope of application	1
ZKS - I - 6L type coagulation test mixer	
Host structure	2
Product Introduction	3
Basic parameters	3
Packing List	3
Product features	4
Product installation	4
Programming Methods	5
<i>Login and password change</i>	5
<i>Writing programs</i>	5
Runs automatically	6
The beep function to turned on/off	7
Loop function to turn on/off	7
Manually run	8
After sales service	9

Preface

Thank you for purchasing and using the coagulation test mixer product from Sichuan Zhongkeshen Technology Co., Ltd!

Before using this product, please carefully read the corresponding user manual and use it in the correct order. Please also keep it properly for future reference.

Special note:

1. This manual is only for use with the corresponding model of product.
2. Please use it according to the specific product model!
3. Our company reserves the final interpretation right of this manual. Thank you again for choosing our product!

Product Overview

The coagulation test mixer is a commonly used water treatment specialized testing instrument in the field of water treatment. It is used for beaker tests in the water treatment process, simulating the processes of dosing, reaction, flocculation, and precipitation in water treatment operations (tap water purification or sewage treatment). The test results are amplified and applied to production equipment, playing an important guiding role in practical production in water plants, helping water plants to accurately control the dosage of chemicals and reduce production costs.

Main role


1. Compare the coagulation effects of various coagulants;
2. Determine the optimal dosage of coagulant;
3. Optimize mixing conditions;
4. Optimize flocculation conditions;
5. Seeking a reasonable combination of mixing, flocculation, and precipitation.

Scope of application

Higher education institutions, research institutes, water treatment plants, sewage treatment plants, municipal water supply and drainage, environmental protection industry.....

ZKS - I - 6L type coagulation test mixer

Host structure

①:Main control touch screen

②:Dosing tube

③:Temperature sensor

④:Mixing cup

⑤:Bottom light board

⑥:Power switch

⑦:Paddle number

⑧:Dosing shaft

⑨:Sampling valve

⑩:Cooling air outlet

⑪:Power interface

⑫:Power cord

⑬:Fuse core

Product Introduction

This product is developed and launched by Sichuan Zhongkeshen Technology Co., Ltd. The water treatment test device adopts a 7-inch touch color LCD screen design, centralized management operation and display, and sensitive response Visually intuitive and clearer. Stepper motor drive, diversified operation modes, dosing, lifting, G value, GT value and other parameters, large capacity storage, program storage up to 20, built-in instrument fault alarm reminder, alarm parameters are clear at a glance.

Basic parameters

Model	ZKS- I -6L
Host weight	27kg
Product size	85*23*62cm
Mixing cup	1L,φ10cm,H16cm
Speed range	10~1000 r/min ±0.01%
Velocity gradient (G)	10~1000 S ⁻¹
Running time	Each program is divided into 10 segments, with each segment running for 0~9 hours
Temperature range	0℃ ~ 50℃ ±1℃
Power	150W
Voltage	AC 220V ±5%

Packing List

Number	Name	quantity	Number	Name	quantity
1	Host	1	5	Mixing cup	6
2	Power cord	1	6	Temperature sensor	1
3	Dosing tube	6	7	Manual	1
4	Fuse core	2	8	Warranty card	1

Product Features

1. Calculate the G value based on the rotational speed;
2. Automatically calculate GT value;
3. Conduct experiments automatically according to programming, add medication automatically, and customize the dosing time and precipitation time;
4. The equipment is equipped with high-quality temperature probes for automatic temperature measurement, with a temperature range of 0 °C~50 °C;
5. The base is equipped with LED light source to ensure brightness and clearer observation of flocculation effect;
6. Several mixing shafts can operate synchronously or asynchronously;
7. Touch color LCD screen, sensitive to touch, intelligent and convenient, dynamically displays various parameters, and data is clearer;
8. It can store 20 of programs, with each program running in 10 stages of continuously variable speed;

Product installation

1. This instrument can only use circular beakers.
2. Take out the mixing cup and place it in the groove of the device base;
3. Remove the medication test tube and install it into the test tube clamp (already installed at the factory);
4. Remove the power cord and insert it into the power interface. After connecting to the power supply, press the "-" direction of the power switch to turn it on. Press the "o" direction to turn off the device. The power supply voltage is AC 220V ± 5%, and the power consumption is under 150W. The power supply must have a reliable grounding wire and two-phase power supply is never allowed. There is a fuse in the power socket on the back of the chassis. If it cannot be turned on normally, you can check if the fuse is burned out;
5. Remove the temperature sensor and insert it into the sensor hole on the right side of the instrument for direct use.

Programming methods

1. Login and password change

1.1 Login: After waiting for a few seconds to turn on the computer and select the system language, enter the main interface, click the "Parameters" button, enter the password (initial password: 111111) to enter the parameters interface, as shown in Figure 1;

1.2. Change Password: Click on "Change Password" in the parameter interface to change the password.


Figure 1: Account Information Interface

2. Programming

Attention: Please confirm that the temperature display is normal before starting to write the program.

Click on the program number box to select the program number you want to edit. In the corresponding program, you can set a total of 10 programs and fill in the "Speed" and "Time" for each stage. The program can automatically calculate the G value based on the speed.

- StirerQuan: Number of mixing heads controlled by the program
- Gd: Calculation coefficient of G value
- Tk, Tb: Temperature coefficient

Fill in the program and click "save" to use it during runtime.

NO.	1	2	3	4	5	6
ProgramNo.	0	0	0	0	0	0
PhaseNo.	0	0	0	0	0	0
Speed	0	0	0	0	0	0
G_Value	0	0	0	0	0	0
GT_Value	0	0	0	0	0	0

AutoOperatio None
 RunTime(S) 0
 Temp 0.0

Figure 2: Auto interface

NO.:Paddle number;

ProgramNO.: The program number of the stirring paddle operation;

PhaseNO.: At which stage of the program the stirring paddle is running;

Speed: The current speed of the stirring paddle;

G value, GT value: The current G value and GT value parameters;

AutoOperatio: Running the current state of the process;

RunTime: How long the program has been running;

Temp: Real-time temperature detected by the temperature sensor.

Runs automatically

1. After setting the operation parameters, the program can start to run synchronously/asynchronously. Click the "Auto" button on the main page to enter the automatic interface, as shown in Figure 2;
2. Click the box corresponding to the program number of the stirring paddle and set the preset program number for the corresponding stirring paddle (if you need to synchronize operation, all the numbers can be set the same);
3. Click parameters to enter the setting interface, as shown in Figure 3. "Dosing time" is the number of seconds in which the program runs to automatically dosing, which can be directly modified by clicking the data box (the maximum of four dosing times are preset, if you need to dosing more times, please use manual dosing, when you need to do so, click "dosing");
4. Click on the "chem_prec_time" box to set the chem precipitate time and return to the automatic interface for the next step;


Figure 3: Parameters interface

- When the automatic interface program Settings are correct, you can click the "Start" button, the program will enter the automatic operation state. When the program is completed, the stirring paddle will automatically rise to the initial position for the precipitation stage. After the precipitation is completed, the device will sound a buzzer, and "Confirm completion" will pop up in the middle of the interface. Click to confirm the buzzer to close;
- During operation, if you need to stop the program, click the "Stop" button to stop the agitator. Wait for a few seconds and the agitator will automatically rise back to its initial position;

The beep function to turned on/off

After the test is completed, a buzzer will sound. Click the buzzer Prompt button to enable or disable the buzzer. If "Buzzer on" is displayed, the buzzer on state is enabled; otherwise, the buzzer is off state.

Loop function to turn on/off

After this function is enabled, the pre-set program will run in a loop for manufacturer testing. If necessary, click the "cyclic" button to turn on and off the cyclic function. If "cyclic on" is displayed, the function will be turned on; otherwise, the function will be turned off, and it will be considered as off by default.

When the loop function is turned on, clicking "stop" on the "automatic monitoring" page will not stop the program and will automatically start the program from scratch.


Figure 4 Manual interface

Manually run

This program is used to test whether the various functions of the equipment are normal, and to verify whether the device can be used normally. This program is used to test whether the various functions of the equipment are normal, and to verify whether the device can be used normally. The left side controls the stirring paddle and the right side controls the dosing shaft.

1. "ORG": Both the stirring paddle and Dosing shaft have default zero positions. After clicking the "ORG" button, the stirring paddle and Dosing shaft will automatically reset;
2. "UP" and "DOWN" control the rise and fall of the mixing paddle, long press to make the mixing head rise/fall;
3. "LifPosition" and "LowerPosition" control the stirring paddle one key to the initial position and one key to the cutoff position;
4. "RotationCW" and "RotationCCW" control the rotation of the dosing axis, and long press to realize the rotation of the dosing axis;
5. "InitialPosition" and "DosingPosition" control the drug axis to return to the preset initial position and pouring position;
6. The "Stir" button can detect whether the corresponding stirring paddle of "NO." is running normally, and the rotation of the stirring paddle can be controlled by long press;
7. "Speed" : control the detection speed of mixing paddle corresponding to "NO.";
8. "Parameters" : Do not modify this parameter without the help of our engineers to avoid abnormal operation of the device.

After-sale service

The company has established a unified standard for after-sales service guidelines, and after-sales engineers have undergone systematic and professional training to provide users with the most professional and caring services.

1. Belonging to product quality issues, replacement within two months and free repair within one year. Material costs will be charged for replacing parts within three years, and maintenance fees will be waived.
2. When there is a quality issue with the product, please contact our company by phone in a timely manner and provide detailed information about the product. Our company will serve you promptly.
3. When there is a quality problem with the product, it may be repaired by returning it to the factory. Please keep the outer packaging of the product.

官方网站
Official website

微信公众号
Wechat public account

网址\Website: www.sczks.com
售后电话: 19138459975
Tel: +8619138459975
服务电话: 028-85541992
Service Calls: +86-28-85541992