MIPARTM

Installation Instructions For Windows

安装包下载网址:	
https://mipar-installer-tokyo.s3.amazona	aws.com/MIPAR v5.1.0.0 Installer Win.zip
一、安装	Setup - MIPAR v4.3.1 - X
1.双击安装包	License Agreement
2. 点击右侧的"I accept the	Pease read die rodowing important mormatori berore continuing.
agreement", 点击"Next"	Please read the following License Agreement. You must accept the terms of this agreement before continuing with the installation.
	IMPORTANTI READ CAREPULT: Inis End-User License Agreement (EULA) is a legal agreement between you (either an individual or a single entity) and MIPAR Software LLC (MIPAR) for the use of certain Software purchased by you, as licensee, from MIPAR, pursuant to an applicable Order Form. An Order Form means an online, electronic or hardcopy communication by or on behalf of you to MIPAR, which contains details of you and the relevant order, including your name, the Software (including whether under a Single Stet License or a Network License) and the number of users (if applicable under a Network License), and for the purpose of downloading, ordening or confirming orders of the Software. Order Forms include the information webpage completed by a licensee on the MIPAR website, or any customer order form (as provided by MIPAR).
	I accept the agreement I do not accept the agreement
	Next > Cancel
3 默认安装路径(切勿修改安装路	Setup - MIPAR v43.1 − □ ×
7. 「古击"Next"	Select Destination Location Where should MIPAR be installed?
	Setup will install MIPAR into the following folder.
	To continue, click Next. If you would like to select a different folder, click Browse.
	C:/Program Files/MIPAR Browse
	At least 2.24 GB of free disk space is required.
	< Back Next > Cancel
	(*) Setup - MIPAR v4.3.1 - X
4.选中"Create a desktop shortout"创	Select Additional Tasks Which additional tasks should be performed?
建桌面快捷方式,点击"Next"	Select the additional tasks you would like Setup to perform while installing MIPAR, then click Next.
	Additional shortcuts:

< Back Next > Cancel

5.点击"Intall"安装,安装需要等待一段时间,点击"Finish"完成安装程序



二、激活

1.双击安装的 **MIPAR** 软件, 试用版请点击 "**Trial**"

2.点击 Activate,选择从邮箱下载的 C2V 密钥文件,或者输入体验版许可证密钥 3.点击 Finish,完成激活程序

Activate						
	T rial	Local	Network	U pdate		
0	Enter Pro	duct Key	Activa	ie.		
	OR					
	Choose License Key Choose the V2C file downloaded from your trial license email					
0	Finish Act	ivation				

MIPARTM

--Spotlight 安装教程--

一、安装

注: MIPAR 软件要求版本: V4.5.0 版本以上。如有需要,请联系科普乐工程师。

1. 安装包下载地址: https://www.keple.cn/MIPAR.html(科普乐官网下载)

2. " 双击 "打开安装包,默认 C 盘。点击	🕑 Setup - MIPAR Spotlight — 🗆 🗙
"Next"下一步。	Select Destination Location Where should MIPAR Spotlight be installed?
	Setup will install MIPAR Spotlight into the following folder.
	To continue, click Next. If you would like to select a different folder, click Browse.
	G1Program Files/MIPAR Browse
	At least 7.57 GB of free disk space is required.
	Next > Cancel
3. 点击"Instal1"安装。	
😕 Setup - MIPAR Spotlight — 🗆 🗙	😢 Setup - MIPAR Spotlight — 🗆 🗙
Ready to Install Setup is now ready to begin installing MIPAR Spotlight on your computer.	Installing Please wait while Setup installs MIPAR Spotlight on your computer.
Click Install to continue with the installation, or click Back if you want to review or change any settings.	Extracting files C:\Program Files\MIPAR\plugins\Spotlight\spotlight_accurate.pt
C:\Program Files\MIPAR	
*	
4	
< Back Install Cancel	Carriel
	Guircu

4. 点击"Finish"完成。



二、激活

注: 要求 MIPAR 软件 V4.5.0 及以上版本; 下载提供的 V2C 许可证; 安装好 Spotlight。

1. 双击安装的 MIPAR 软件。

2. 打开任意模块,点击"Help"->"Activate License"->"Update"。

3. 单击 "Activate",选择邮箱下载的 "Spotlight 的 V2C" 文件。

4. 重新启动 MIPAR 软件, Spotlight 和 Snap 已经安装激活完毕。加载图形就可 正常使用。



注: 一台电脑不能试用第二次, 如有问题, 请联系客服: 024-31014540